



Transverse Maxillary Distraction


Rotterdam Palatal Distractor



Rotterdam Palatal Distractor
1:1 scale



Rotterdam Palatal Distractor, start position



During distraction period



Rotterdam Palatal Distractor

Distractors	Item Number
9 mm initial width (19 mm distraction length)	51-555-09-09 Top Seller
13 mm initial width (19 mm distraction length)	51-555-13-09

Distraction length/turn see distraction diagram


Screws are not required!

Patient screwdrivers

Hockey stick-like	51-555-90-07
Straight	51-555-85-07
Angled	51-555-95-07

Distraction Diagram

showing the width of a 9-mm Rotterdam Palatal Distractor in relation to active distraction time




For additional product information, please refer to the brochure "Rotterdam Palatal Distractor".

Rapid Palatal Expander



Rapid Palatal Expanders (RPEs)



RPE Distractor



51-565-90-07



2 maxDrive® Drill-Free screws
2.0 x 7 mm



Additionally available:
Flexible activating wrench
51-565-95-07

Rapid Palatal Expander (RPE)

Distractors	Item Number	sterileR
① 9 mm distraction length	51-565-09-09*	51-565-09-71**
② 18 mm distraction length TopSeller	51-565-18-09*	51-565-18-71**
③ 27 mm distraction length	51-565-27-09*	51-565-27-71**
④ 36 mm distraction length	51-565-36-09*	51-565-36-71**

Recommended distraction length
1-2 color codes = 0.33 - 0.66 mm/day
(one complete turn = 1.0 mm)

Measuring Templates


Size I	51-564-09-09
Size II	51-564-18-09
Size III	51-564-27-09
Size IV	51-564-36-09

Recommended screws

Drill-Free: 2.0 x 7 mm

Patient screwdriver

Activating wrench	51-565-90-07
Flexible activating wrench	51-565-95-07



1:1 scale

Δ max. extension

* distractor including activating wrench 51-565-90-07 and 2 maxDrive® Drill-Free screws 2.0 x 7 mm
** sterile packed distractor, including activation wrench 51-565-90-07 and 2 maxDrive® Drill-Free screws 2.0 x 7 mm